

ISO80369-6非対応製品と代替品一覧 (2019年11月12日作成)

No.	ISO80369-6非対応製品			代替品		
	JANコード	製品コード	サイズ、型式等	JANコード	製品コード	サイズ、型式等
	販売名称：ディスプレイザブル神経ブロック針 (一般的名称：麻酔用滅菌済み穿刺針)			販売名称：ディスプレイザブル神経ブロック針 MM (一般的名称：麻酔用滅菌済み穿刺針)		
1	4526737270753	27112010	20Gx40mm	-	-	-
2	4526737270777	27112030	20Gx60mm	4526737275727	27112036	20Gx60mm
3	4526737270784	27112040	20Gx80mm	4526737275734	27112046	20Gx80mm
4	4526737270791	27112050	20Gx100mm	4526737275741	27112056	20Gx100mm
5	4526737270807	27112060	20Gx120mm	4526737275758	27112066	20Gx120mm
6	4526737270814	27112070	20Gx140mm	4526737275765	27112076	20Gx140mm
7	4526737270821	27112110	21Gx40mm	-	-	-
8	4526737270838	27112120	21Gx60mm	4526737275772	27112126	21Gx60mm
9	4526737270845	27112130	21Gx80mm	4526737275789	27112136	21Gx80mm
10	4526737270852	27112140	21Gx100mm	4526737275796	27112146	21Gx100mm
11	4526737270869	27112150	21Gx120mm	4526737275802	27112156	21Gx120mm
12	4526737270876	27112160	21Gx140mm	4526737275819	27112166	21Gx140mm
13	4526737270883	27112200	22Gx40mm	4526737275826	27112206	22Gx40mm
14	4526737270906	27112230	22Gx60mm	4526737275840	27112236	22Gx60mm
15	4526737270920	27112250	22Gx80mm	4526737275864	27112256	22Gx80mm
16	4526737270937	27112260	22Gx100mm	4526737275871	27112266	22Gx100mm
17	4526737270944	27112270	22Gx120mm	4526737275888	27112276	22Gx120mm
18	4526737270951	27112280	22Gx140mm	4526737275895	27112286	22Gx140mm
19	4526737270975	27112310	23Gx40mm	4526737275925	27112316	23Gx40mm
20	4526737270982	27112320	23Gx60mm	4526737275932	27112326	23Gx60mm
21	4526737270999	27112330	23Gx80mm	4526737275949	27112336	23Gx80mm
22	4526737271002	27112340	23Gx100mm	4526737275956	27112346	23Gx100mm
23	4526737271019	27112350	23Gx120mm	4526737275963	27112356	23Gx120mm
24	4526737271026	27112360	23Gx140mm	4526737275970	27112366	23Gx140mm
25	4526737271057	27112420	24Gx40mm	4526737275994	27112426	24Gx40mm
26	4526737271040	27112410	24Gx60mm	4526737275987	27112416	24Gx60mm
27	4526737271064	27112510	25Gx40mm	4526737276007	27112516	25Gx40mm
28	4526737271071	27112520	25Gx60mm	4526737276014	27112526	25Gx60mm
29	4526737271088	27112530	25Gx80mm	-	-	-
30	4526737271095	27112540	25Gx100mm	4526737276021	27112546	25Gx100mm
31	4526737271101	27112550	25Gx120mm	-	-	-
	販売名称：スパイナル針 M型 (一般的名称：麻酔脊髄用針)			販売名称：スパイナル針 P型 (一般的名称：麻酔脊髄用針)		
32	4526737270180	27031810	18Gx3”(76mm)	4526737275192	27011876	18Gx70mm
33	4526737270197	27031820	18Gx3・1/2”(90mm)	4526737275208	27011896	18Gx90mm
34	4526737270203	27031910	19Gx3”(76mm)	4526737275215	27011976	19Gx70mm
35	4526737270210	27031920	19Gx3・1/2”(90mm)	4526737275222	27011996	19Gx90mm
36	4526737270227	27032010	20Gx3”(76mm)	4526737275239	27012076	20Gx70mm
37	4526737270234	27032020	20Gx3・1/2”(90mm)	4526737275246	27012096	20Gx90mm
38	4526737270258	27032110	21Gx3”(76mm)	4526737275253	27012176	21Gx70mm
39	4526737270265	27032120	21Gx3・1/2”(90mm)	4526737275260	27012196	21Gx90mm
40	4526737270302	27032210	22Gx3”(76mm)	4526737275277	27012276	22Gx70mm
41	4526737270319	27032220	22Gx3・1/2”(90mm)	4526737275284	27012296	22Gx90mm
42	4526737270357	27032310	23Gx3”(76mm)	4526737275314	27012376	23Gx70mm
43	4526737270364	27032320	23Gx3・1/2”(90mm)	4526737275321	27012396	23Gx90mm
44	4526737270456	27032510	25Gx3”(76mm)	4526737275390	27012576	25Gx70mm
45	4526737270463	27032520	25Gx3・1/2”(90mm)	4526737275406	27012596	25Gx90mm
46	4526737270494	27032610	26Gx3”(76mm)	-	-	-
47	4526737270500	27032620	26Gx3・1/2”(90mm)	-	-	-

ISO80369-6非対応製品の針基は金属ですが、代替品の針基はプラスチックです。

ISO80369-6非対応製品の針基は金属ですが、代替品の針基はプラスチックです。

また、長さ76mmのISO80369-6非対応製品の代替品の長さは70mmとなります。

※ISO80369-6非対応製品である『ディスプレイザブル神経ブロック針』『スパイナル針 M型』はいずれも2020年2月末をもちまして販売を終了させていただきます。